

# Cover-URO®

*Disposable sterile sheaths for flexible cystoscopes*

- Sterile, practical and ready to use
- Enhances safety & productivity
- Less risk of infection and adverse events



- Preserves efficiency of the flexible cystoscope
- Reduces service expenses and losses of time
- Effective barrier against viruses and bacteriums

**C+VER**

# Cover-URO®

*Disposable sterile sheaths for flexible cystoscopes*

**Cover-URO®**

is a sterile sheath used as a protection of flexible cystoscopes during cystoscopy procedures. Compatible with the most commonly used instruments in urologic field, it assures sterile and fast examinations.



## Cystoscope Model

STORZ Mod. 11272C-11272VP-11272VN

OLYMPUS Mod. CYF-5/CYF-5A/CYF-5R

WOLF Mod. 7305.001/7305.006

PENTAX Mod. ECY-1575K/ECY-1570

## Product Code

GCSSTO01

GCSOLY01

GCSWOL01

GCSPEN01

## Packaging:

A dispenser contains 10 individually packaged sheaths.

A box contains 5 dispensers.

Produced by:

**COVER**

C.so V. Emanuele II, 12 – 10123 Torino (Italia)

Tel 0182 077213 – Fax 0182 989819

email: [info@cover-srl.com](mailto:info@cover-srl.com)

web site: [www.cover-srl.com](http://www.cover-srl.com)



0476  
UNI CEI EN ISO  
13485:2012

Distributed by: