

DISPOSABLE COVER FOR HYSTEROSCOPY


Gynko[®]

Ø 3mm

with 5 or 7 fr. CHANNEL



- Utilises a patented single use 100% sterile Gynko® Sheath
- Minimal trauma & speedy procedure for patient - 2.9mm Ø optic
- Multiple patients procedures in one session
- Perform minor operative procedures using 40cm flexi scissors / forceps
- No requirement for external decontamination services - significant cost savings
- No risk of damage / wear & tear of hysteroscopy equipment by outsourcing to clean



GFRGYN300

GYNKO® for rigid optics, 2,9mm, 300mm, 0° rigid lens with Storz-like connection and camera cover, 5 fr. channel **Box: 10 pieces**

GFRGYN300W

Without camera cover **Box: 10 pieces**

GFRGYN330

GYNKO® for rigid optics, 2,9mm, 300mm, 30° rigid lens with Storz-like connection and camera cover, 5 fr. channel **Box: 10 pieces**

GFRGYN330W

Without camera cover **Box: 10 pieces**

GFRGYN337

GYNKO® for rigid optics, 2,9mm, 300mm, 30° rigid lens with Storz-like connection and camera cover, 7 fr. channel **Box: 10 pieces**

GFRGYN337W

Without camera cover **Box: 10 pieces**

Gynko®

Distributed by:

